DURONIC
Creative Consumer Products

OX02
Fingertip Pulse Oximeter

PLEASE READ ALL THE INSTRUCTIONS BEFORE USE

duronic.com duronic @duronicmedia Duronic @duronicmedia

Contents

English Manual	**3**
Manuel d’utilisation en français	**19**
Gebrauchsanleitung auf Deutsch	**32**
Manual de Instrucciones en Español	**45**
Manuale D’istruzioni In Italiano	**58**
Instrukcja obsługi. Polski	**71**
Türkçe Kullanım Kılavuzu	**84**
Gebruikershandleiding Nederlands	**97**
Svensk manual	**110**

English Manual

DURONIC

Introduction


Thank you for purchasing the Duronic OX02 Fingertip Pulse Oximeter. This Manual is written and compiled in accordance with Medical Device Regulation (MDR) 2017/745 of the European Parliament and of the Council. In case of modifications and software upgrades, the information contained in this document is subject to change without notice.

Please read this User Manual carefully before using this product and retain for future reference.

Box Contents

- 1 x Oximeter
- 2 x AAA batteries
- 1 x Lanyard
- User Manual

Digital Display

1. Pulse Oxygen Saturation (SpO_2)
2. Pulse Rate (PRbpm)
3. Perfusion Index (PI%) + Respiration Rate (RR/min)
4. Pulse Volume Wave
5. Battery Indicator
6. Pulse Bar

Note: The blood perfusion index (PI%), and respiratory rate (RR) only show status values and are not used for clinical reference.

DURONIC

Safety Information

PLEASE READ THIS MANUAL CAREFULLY AND RETAIN FOR FUTURE REFERENCE

Instructions for Safe Use

- Check the main unit and all accessories periodically to make sure that there is no visible damage that may affect patient's safety. If used regularly, it is recommended that the device should be inspected once a week at least. When there is obvious damage, stop using the monitor.
- You must not attempt to repair or disassemble the oximeter yourself. If the device requires repair, please contact the Duronic Customer Care Team.
- The oximeter cannot be used together with devices not specified in User Manual. The only accessories (if any) that can be used with this device are those recommended by manufacturer.
- This oximeter is calibrated before leaving factory.

⚠ **Warning:**

- Explosive hazard---DO NOT use the oximeter in environment with inflammable gas, for example anaesthetic mixture with air/oxygen or where there is nitrous oxide. These are highly flammable.
- DO NOT use the oximeter in conjunction with MRI and CT equipment.
- This device is not suitable for anyone who is allergic to ABS or silicone.
- You may feel discomfort or pain if using the device continuously. It is recommended that the sensor should not be applied to the same finger for over 2 hours at a time.
- The infrared light emitted from the device is harmful to the eyes. Do not stare at the light.

Safety Information

Attention

- Please check the contents before use to make sure the device and accessories are included in accordance with the list in this manual. If you are missing any parts, please contact Duronic directly for assistance.
- Keep the oximeter away from dust, vibrations, corrosive substances, explosive materials, high temperatures and moisture.
- Please keep the product out of reach from children and pets. It is not a toy.
- Enamel, Gel or Acrylic fingernail polish or other fingernail applications may distort and/or produce inaccurate readings. Please remove them before use.
- If the oximeter gets wet, stop using it immediately.
- When moved from a cold environment to a warm or humid environment, please do not use the device straight away.
- Do not expose to high temperature or clean with high pressure steam as this will damage the oximeter.
- Do not submerge the oximeter in water or any other liquid. When it needs cleaning, please wipe its surface with medical alcohol with a soft cloth. Do not spray any liquid on the device directly.
- Do not use this device on young infants or newborn babies.
- This oximeter is suitable for children above four years old and adults (weight should be between 15kg to 110kg).
- The device may not work for all patients. If you are unable to achieve stable readings, discontinue use.
- The update period of data is less than 5 seconds, which is changeable according to an individual’s pulse rate.
- The oximeter does not have low-battery alarm function, it only shows the low-battery symbol on the digital display. Please change the battery when the energy has run out.
- Batteries must be removed if the device is going to be stored for more than one month, otherwise the batteries may expire and leak.
- The disposal of the oximeter, its accessories and/or packaging (including battery, plastic bags, foams and paper boxes) should be done so responsibly following the local council/district refuse and recycling regulations.

About the OX02

What Does it Do

The OX02 oximeter measures two things: oxygen saturation and pulse rate.

- Oxygen saturation is the percentage of oxyhaemoglobin (HbO_2) in blood. HbO_2 is a combination of haemoglobin (Hb) and oxygen molecules (O_2), that is present in all human blood and is an important physiological factor involved in respiration and circulation. The oxygen saturation of arterial blood in the average human body is around 98%. Measuring oxygen saturation is an important indicator of the amount of oxygen in the body. If there is not enough oxygen it can affect the normal function of the body in many ways and result in a condition called hypoxia.

- In general, the normal levels of oxygen saturation should not be lower than 94%; if a measurement of 94% or less has been obtained, it means an insufficient supply of oxygen has been detected and proper medical advice should be sought from a qualified doctor.

- Pulse rate is the number of pulse beats per minute. In general, the pulse rate of every person should be between 60-100 beats per minute. If you are concerned about your pulse being too fast or slow, consult a qualified doctor for advice.

Note:

- Frequent heavy exercise can affect your normal pulse rate to register lower than the standard person. Tell your doctor if seeking medical assistance or advice.

Note:

- If unsure on how to use the device correctly for a medical reason, talk to your doctor before using a home device like the OX02 so that you understand how to interpret the results.

About the OX02

How Does it Work?

The Pulse Oximeter uses light to work out oxygen saturation. Light is emitted from one side of the device which goes across the pulse oximeter probe and reaches the light detector on the other side. When a finger is placed in the oximeter between the source of light and the detector, the light has to try and pass through the finger and through the blood cells. By measuring how much light reaches the light detector on the other side, the pulse oximeter knows how much light has been absorbed by the haemoglobin (Hb) in the blood. The more Hb in the finger there is, the more light absorbed.

Application of the OX02 Oximeter

- The OX02 Pulse Oximeter is a non-invasive device intended for the spot-check of oxygen saturation of arterial haemoglobin (SpO_2) and the pulse rate of adult and paediatric patients. With a clear digital display, this device is easy to use and comes with a lanyard to carry it with you should you wish.

- The oximeter can be used in a wide range of situations, for instance, family care at home, in hospitals, for social medical care organisations and for sports. While the OX02 can be used before and after sports, it is not recommended to be used during sports activities.

Warning:

This device is not suitable for use in patients requiring continuous supervision.

Operating Guide

Installing the Batteries

- Step 1: Slide the cover to unlock, then lift the cover off.
- Step 2: Insert the two AAA size batteries in the correct position.
- Step 3: Replace the cover ensuring it is locked in place.

Warning:

Take care when you insert the batteries as improper insertion may damage the device.

Attaching the Lanyard

- Step 1: Insert the loop of the lanyard into one of the holes and push the loop in until it comes out of the other hole.
- Step 2: Thread the lanyard strap through the loop and pull tight to secure in place.

Rotating the Display

Press the button once to turn on. Once on, press the button again to turn the display 90°. Keep pressing the button to rotate further. There are 4 positions available.

Parameter Setup

When the readings fall below the parameters, the oximeter will continuously beep. Press and hold the button to enter the alert menu. Short-press the button to scroll through the options or long-press to select. Here you can change the settings;

- to turn on / off the alert indicator select ‘Beep’, then long-press the button until it reads ‘on’ / ‘off’ next to it.

- to change the warning parameters of SpO_2 and PR select the measurement, then long-press the button. Select '+/-' and long-press to change the long-press to ('+') increase or ('-') decrease the value.
- to restore factory settings, select 'Restore', then long-press the button until it reads 'on' next to it.

How to Take a Measurement

Only use the oximeter when you are sitting or standing comfortably. You will not get an accurate reading if you are under exertion or have been doing any vigorous activity before or during use of the device.

1. Insert the two batteries properly and then replace the cover.
2. Open the clip by pinching it open.
3. Place your finger nail-side up into the device resting it on the pads inside. Make sure your finger is fully inserted, ensuring the end of the finger touches the end of the space inside and that the finger's first digit covers the sensor.
4. Close the clip around the finger
5. Press the button on front panel once.
6. Do not shake the finger. Stay still and relaxed during the process.
7. The oximeter will instantly give you a reading, providing you with your current statistics on the display screen.

Maintenance and Cleaning

- Change the batteries when the battery symbol flashes on the screen.
- To clean the device, wipe with damp (well wrung) cloth. Then leave to air dry or dry with a clean dry cloth.
- Medical alcohol can be used in place of water if you wish to disinfect the product after use and prevent cross infection. Please clean between each user.
- If you plan to not use the oximeter for a long period of time, remove the batteries before storing.

⚠ **Warning:**

High-pressure, high-temperature or gas sterilization cannot be used on the device.

Do not immerse the device in water or any other liquid.

It is recommended that the device should be kept in a dry environment. Humidity may reduce the life of the device, or even damage it.

Troubleshooting

Problem	Possible Reason	Solution
The SpO_2 and Pulse Rate is not displaying normally	a. The finger is not properly positioned. b. The patient's SpO_2 is too low to be detected.	a. Place the finger properly and try again. b. Try again; If you are sure that the device is working well, seek medical help.
The SpO_2 and Pulse Rate are not displayed stably	a. The finger is not placed inside deep enough. b. The finger is shaking or the patient is moving.	a. Place the finger properly and try again. b. Calm the patient and then try again.
The device cannot be turned on	a. Low battery or no battery. b. The batteries are not inserted properly. c. The malfunction of the device.	a. Change batteries. b. Check batteries are installed correctly. c. Contact Duronic.
The display has turned off suddenly	a. The device will power off automatically when there is no signal within 8 seconds. b. The batteries are drained.	a. Normal. b. Change batteries.

Technical Specifications

Display Information	**Display Mode (Digital TFT)**
Pulse Oxygen Saturation (SpO_2)	Yellow
Pulse Rate (PRbpm)	Green
Perfusion Index (PI%)	(1 decimal place) - Blue
Respiration Rate (RR/min)	Blue
Pulse Volume Wave	Red
Pulse Bar	Blue
SpO_2 Parameter Specification	
Measuring range	70%-100% (Resolution: 1%)
Accuracy	80%-100%; ±2%, 70%-79%; ±3%
Pulse Parameter Specification	
Measuring range	30bpm-254bpm (Resolution: 1 bpm)
Accuracy	≤100BPM; ±1BPM, >100BPM; ±2BPM
Alert conditions (Factory Settings)	
SpO_2	Less than 94%
PR	Less than 50bpm or more than 130bpm
Allergens	
ABS, Silicone (Latex Free) - If a reaction occurs discontinue use immediately.	
Battery Requirement	
2 X 1.5V (AAA size) alkaline battery	
Battery Useful Life	
Two batteries can work continuously for 6 hours (battery life will vary)	
Power off	
The Oximeter will automatically turn off when no finger is inserted within 8 seconds.	
Optical Sensor	
Red light (wavelength: 660nm) 1.8mW	Infrared (wavelength: 940nm) 2.0mW
Dimensions and Weight	
Dimensions	59 (L) X 32.5 (W) X 29.3 (H) mm
Weight	About 28g (without batteries)

Key of Symbols

	Type BF applied part		Keep away from sunlight
	Date of manufacture	SN	Serial number
	Manufacturer		Separate collection for electrical and electronic equipment
	Use-by date	IP22	Protected against solid foreign objects greater than 12.5mm in diameter and dripping water when tilted up to 15º
	Caution		
MD	Medical device		This device has no alarm function
CE 0123	CE mark		Not made with natural rubber latex
EC REP	Authorized representative in the EU		Keep dry
	Mobius Loop		Refer to instruction manual
RoHS	Restriction of Hazardous Substances		Green Dot

DURONIC

Device Specifications

Classification:
Class IIa (MDR) 2017/745
EMC complies with IEC60601-1-2 standard
Anti-electric Shock Type: Internally powered equipment
Anti-electric Shock Degree: Type BF equipment
Overvoltage category classification: Class I
Pollution Degree 2

Environment Requirements

Operation Temperature:	5°C ~ 40°C
Storage Temperature:	-10°C ~ 40°C
Ambient Humidity:	15%-80% RH, no condensation in operation 10%-95% RH, no condensation in storage
Atmospheric Pressure:	70kPa ~ 106 kPa, in operation 70kPa ~ 106 kPa, in storage

Xuzhou Yongkang Electronic Science Technology Co., Ltd
1st & 2nd Floor, 6#01, 6#02 | No.6 Building 1st Phase Economic Development Manufacturing Zone | LIANDO U Valley | No.6 Leye Road | Xuzhou ETDZ | 221000 Xuzhou | People's Republic of China

EC REP **Prolinx GmbH**
Brehmstr. 56 | 40239 Düsseldorf | Germany

MPN: YK-80B

DURONIC **UK Distributor: Shine-Mart Ltd**
1 Spilsby Road | Romford | RM3 8SB | UK

EU Distributor: Shinemart GmbH
Brandenburger Str. 40 | 40880 Ratingen | DE

Appendix

Table 1: Guidance and manufacture's declaration – electromagnetic emissions-for all EQUIPMENT and SYSTEMS

Guidance and manufacture's declaration – electromagnetic emission		
The Duronic OX02 is intended for use in the electromagnetic environment specified below. The customer or the user of the OX02 should assure that it is used in such an environment.		
Emission test	**Compliance**	**Electromagnetic environment – guidance**
RF emissions CISPR 11	Group 1	The OX02 uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emission CISPR 11	Class B	The OX02 is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	N/A	N/A
Voltage fluctuations / flicker emissions IEC 61000-3-3	N/A	N/A

Appendix

Table 2: Guidance and manufacturer's declaration – electromagnetic immunity – for all EQUIPMENT and SYSTEMS

Guidance and manufacturer's declaration – electromagnetic immunity		
The Duronic OX02 is intended for use in the electromagnetic environment specified below. The customer or the user of the OX02 should assure that it is used in such an environment.		
Immunity test	**EN 60601 test level**	**Compliance level**
Electrostatic discharge (ESD) IEC 61000-4-2	± 8 kV contact ± 15 kV air	± 8 kV contact ± 15kV air
Electrostatic transient / burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	N/A
Surge IEC 61000-4-5	± 1 kV differential mode ± 2 kV common mode	N/A
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	< 5 % UT (>95 % dip in UT) for 0.5 cycle 40 % UT (60 % dip in UT) for 5 cycles 70 % UT (30 % dip in UT) for 25 cycles < 5 % UT (>95 % dip in UT) for 5 sec	N/A
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m	30 A/m

Table 3: Guidance and manufacturer´s declaration – electromagnetic immunity – for EQUIPMENT and SYSTEM that are not LIFE-SUPPORTING

Guidance and manufacturer´s declaration – electromagnetic immunity		
The Duronic OX02 is intended for use in the electromagnetic environment specified below. The customer or the user of the OX02 should assure that it is used in such an environment.		
Immunity test	**EN 60601 test level**	**Compliance level**
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	N/A

Appendix

Table 4: Test specifications for enclosure port immunity to RF wireless communications equipment

Test frequency (MHz)	Band[a] (MHz)	Service[a]	Modulation[b]	Max power (W)	Distance (m)	Immunity test level
385	380 -390	TETRA 400	Pulse modulation[b]	1,8	0.3	27
450	430 - 470	GMRS 460, FRS 460	FM[c] ± 5 kHz deviation 1 kHz sine	2	0.3	28
710	704 - 787	LTE Band 13, 17	Pulse modulation[b] 217 Hz	0,2	0.3	9
745						
780						
810	800 - 960	GSM 800/900. TETRA 800, iDEN 820, CDMA 850, LTE Band 5	Pulse modulation[b] 18 Hz	2	0,3	28
870						
930						
1 720	1700 - 1990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3 4, 25; UMTS	Pulse modulation[b] 217 Hz	2	0,3	28
1 845						
1 970						
2 450	2400 - 2570	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation[b] 217 Hz	2	0,3	28
5 240	5100 - 5800	WLAN 802.11 a/n	Pulse modulation[b] 217 Hz	0,2	0.3	9
5 500						
5 785						

NOTE: If necessary to achieve the immunity test level, the distance between the transmitting antenna and the me equipment or me system may be reduced to 1 m. The 1 m test distance is permitted by IEC 61000-4-3.

a) For some services, only the uplink frequencies are included.

b) The carrier shall be modulated using a 50 % duty cycle square wave signal.

c) As an alternative to FM modulation, 50 % pulse modulation at 18 Hz may be used because while it does not represent actual modulation, it would be worst case.

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DURONIC

WEEE

UK **Information for consumers on the disposal of Waste Electrical & Electronic Equipment (WEEE).**

This mark on a product and/or accompanying documents indicates that when it is to be disposed of, it must be treated as Waste Electrical & Electronic Equipment, (WEEE).

Any WEEE marked products must not be mixed with general household waste, but kept separate for the treatment, recovery and recycling of the materials used. For proper treatment, recovery and recycling; please take all WEEE marked waste to your Local Authority Civic waste site, where it will be accepted free of charge.

If customers dispose of Waste Electrical & Electronic Equipment correctly, they will be helping to save valuable resources and preventing any potential negative effects upon human health and the environment, of any hazardous materials that the waste may contain.

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